


FLOOME™ - QUICK START GUIDE

Note: the FLOOME™ app requires an active Internet connection.

1. Download and install the FLOOME™ app. It is available for free from the Apple App Store, Google Play Store and Windows Phone Store.
2. Create a new user account, or login with an existing one (at the first login, you will be presented with a short tutorial and the request to insert your body's info).
3. From the main window, you can tap the “menu” button to access the menu or the “Help” button to open the tutorial (see the following picture)



4. Rotate the silicone mouthpiece and extend the 3.5mm jack from the device.
5. Connect your FLOOME™ to the smartphone, and wait a few seconds to let the app recognize the device. After that, a “Start” message will be displayed on the screen along with 3 blinking arrows. Scroll in the direction shown by the arrows and wait until FLOOME™ is ready to begin the measurement.



6. When “Blow” appears on the screen, start (and keep) blowing into the mouthpiece like if you were whistling. A whistle sound is emitted by FLOOME™ while blowing: try to keep it as constant as possible. When your smartphone’s display and flashlight blink or the phone vibrates, stop blowing (the flashlight may not blink on some smartphone models or if the app wasn’t granted access to the camera).
7. Your blood alcohol content (BAC) will be calculated. The measurement progress is indicated by a white progress wheel. Finally, your BAC will be displayed along with the appropriate legal limit.

WHAT TO AVOID

When using FLOOME™, please follow these simple hints in order to avoid inaccurate results or measure failings:

- ✗ Don't cover the hole on the top shell of the device while blowing.
- ✗ Don't blow into the device before the application tells you to do so. If you did it, wait at least one minute before taking a measurement to allow the sensor to clean up.
- ✗ Don't blow into the device if it is not plugged into the smartphone.
- ✗ Don't blow just after drinking, but wait at least 20 minutes. Otherwise, the result won't be accurate.
- ✗ Don't unplug the device from the smartphone while a measurement is still in progress.

TROUBLESHOOTING

1. The measure failed

You may need a little practice to learn how to blow correctly. The first times you use the device, you could either blow too hard or too soft. In both cases, the measurement will fail and this will be notified by the application with an error message: "Test failed". With a tap on the error message, you can obtain a more detailed description of the error, along with a suggestion for fixing it.

2. The measure failed so why should I wait 1 min. before taking another test?

Every time you blow, the sensor inside FLOOME™ needs some time to recover, even if the measurement fails. A measurement performed while the sensor has not fully recovered may lead to an inaccurate result. For this reason, you always need to wait at least one minute between two consecutive tests, regardless the outcome of the measure.

FAQ (Frequently Asked Questions)

The FAQs are available on the FLOOME™ website, at the address <https://www.floome.com/faq>.